



Oxford Performance Materials

30 South Satellite Road • South Windsor, CT 06074 USA • T: +1.860.698.9300 • F: +1.860.698.9978 • www.oxfordpm.com



High Performance Thermoplastic OXPEKK® Materials

OXPEKK® BioFlex™ Laminates for Implant Applications

OXPEKK® BioFlex™ smart materials are composed of multiple layers of continuous carbon fibers and thermoplastic OXPEKK® resin. The fibers are orientated within the part, so that stiffness is defined by fiber direction. The isotropic nature of metallic and conventional polymeric materials are overcome, thus opening a new era of dynamic stabilization and motion preservation.



OXPEKK® BioFlex™ Cross-Ply

To date orthopedic implant design engineers have relied upon complex mechanical assemblies and articulating surfaces in order to preserve motion. OXPEKK® BioFlex™ introduces a new and dynamic tool for the orthopedic OEM's via the combination of smart

materials and established design techniques.

OXPEKK® BioFlex™ can be supplied in oriented laminated plates of thickness 1 – 10mm up to 60cm square, semi-finished thermoformed 3D parts waterjet cut and hot-formed from oriented plates and solid unidirectional rods. Common orientations include Cross-Ply (0/90), Unidirectional, Quasi-Isotropic (-45/0/45/90), Torsional-Ply (-45/45). Parts can be readily machined on conventional equipment such as is used for other carbon fiber composites.



OXPEKK® BioFlex™ Torsion-Ply

OXPEKK® BioFlex™ is available for long-term implants with an executed supply agreement.

Mechanical Properties			
0° Tensile		90° Tensile	
Strength MPa/ksi	2506/364	Strength MPa/ksi	57/8.3
Modulus GPa/Msi	146/21.1	Modulus GPa/Msi	10.7/1.6
0° Compressive		90° Compressive	
Strength MPa/ksi	1615/234	Strength MPa/ksi	238/34.5
Modulus GPa/Msi	123.2/17.9	Modulus GPa/Msi	10.7/1.55
Open Hole		Open Hole	
Tensile MPa/ksi	398/65.6	Compressive MPa/ksi	335/50.5

NOTICE TO USERS: To the best of our knowledge, the information contained in this publication is accurate, however we do not assume any liability whatsoever for the accuracy and completeness of such information. The analysis techniques included in this provisional data sheet are often simplifications and, therefore, approximate in nature. Any determination of the suitability of a particular material for any use contemplated by the user is the sole responsibility of the user. Material data and values included in this provisional data sheet are either based on testing of laboratory test specimens and represent data or were extracted from various published sources. All are believed to be representative. These values are not intended for use in establishing maximum, minimum, or ranges of values for specification purposes. It is the sole responsibility of the users to investigate whether any existing patents are infringed by the use of the materials mentioned in this provisional data sheet.